


1 Log in using a 'Member' access

2 Enter the title of the procedure to be customized in the search engine, then select it

- See the “Searching for a procedure using keywords” Tool if necessary


3 Click on the “Customize” button located on the right of the procedure header

Insertion and removal of a suprapubic urinary catheter



4 Select the “Tools” tab

5 Click on “Choose Files” in the “Add a new file” box



6 Select the file(s) to upload in the ‘Tools’ section of the platform, then click “Open”

- Check that each uploaded file does not exceed the 100 Mo limit
- Renaming the files before uploading them is recommended so that the file names are those that will be displayed

7 Click on “Save” at the bottom left of the page

